SNAP[®] Bile Acids Test Kit

For use only with the IDEXX SNAP[®] Reader Analyzer and the IDEXX SNAPshot Dx[®] Analyzer.

For veterinary use only.

English Version

Overview


The SNAP Bile Acids Test Kit is an enzyme-linked immunosorbent assay (ELISA) for the quantitative measurement of bile acids in dog or cat (domestic) serum. For in vitro use only.

SNAP Bile Acids Test Ranges

Dynamic range: 5–30 µmol/L

Precautions and Warnings

- Run the SNAP calibration device at the beginning of each day (SNAP Reader analyzer only).
- This test cannot be visually interpreted. The IDEXX SNAP Reader analyzer or SNAPshot Dx analyzer must read all tests.
- Use a separate pipette tip and sample tube for each test.
- Do not use components after their expiration dates.
- Do not mix components from different lots of test kits.
- Do not use a device that does not have two blue stripes visible before adding the sample.
- Keep the device in its sealed foil pack until ready to use.
- Keep the SNAP device horizontal during the entire test procedure.
- IMPORTANT: Do not depress the activator until instructed.**
- Insert the device into the analyzer immediately after activation.
- WARNING:** Conjugate – H316/P332+P313/EUH208. Causes mild skin irritation. If skin irritation occurs: Get medical advice/attention. Contains Kathon. May produce an allergic reaction.

Storage

SNAP devices and test reagents must be stored at 2°–7°C (36°–45°F). Bring components to room temperature (18°–27°C/64°–80°F) before running tests. This will take at least 30 minutes, depending upon the temperature in your laboratory.

- Do not heat or microwave the kit components.
- Take out only the kit components you will need.
- Do not leave kit components at room temperature for more than 8 consecutive hours.

Kit Components

- Conjugate solution
- SNAP test devices
- Disposable pipette tips for the 300-µL manual pipettor
- Disposable pipette tips for the 100-µL manual pipettor (for use with the SNAPshot Dx analyzer only)
- Disposable sample tubes
- Workstation

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
Components Required (Not Provided in Kit)

SNAP Reader analyzer users:

- VetTest[®] pipette tips
- 300-µL manual pipettor (provided with the SNAP Reader analyzer)

SNAPshot Dx analyzer users:

The following required components are provided with the SNAPshot Dx analyzer:

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Sample Information

- Do not use plasma or whole blood. Use serum, either fresh, previously frozen or stored at 2°–7°C (36°–45°F) for up to 7 days. To store longer than 7 days, freeze the sample (-20°C/-4°F or below). Remove serum from the red blood cells before freezing. Bring previously refrigerated or frozen samples to room temperature (18°–27°C/64°–80°F) before running assays.
- SNAP bile acids sample size requirements: 100 µL.
- Centrifuge previously frozen or older samples before use.
- Results of this test are not affected by hemolysis or lipemia.

Test Procedure

SNAP Reader analyzer users:

The VetTest analyzer provides on-screen instructions throughout the analysis. For complete instructions, refer to the Running Quantitative SNAP Tests chapter in the SNAP Reader Operator's Manual.

SNAPshot Dx analyzer users:

The SNAPshot Dx analyzer provides on-screen instructions throughout the analysis. For complete instructions, refer to the Using the SNAPshot Dx Analyzer chapter in the SNAPshot Dx Operator's Manual.

Test	Sample	Conjugate
Preprandial	100 µL	300 µL
Postprandial	100 µL	300 µL

Interpreting Test Results

The recommendations below are intended to provide general guidance only. As with any diagnosis, you should use clinical discretion with each patient based on a complete evaluation of the patient, including history, physical examination and complete laboratory profile.

	Postprandial Result			
	<12 µmol/L	12–25 µmol/L	>25 µmol/L	
Preprandial Result	<12 µmol/L	Normal; consistent with sufficient hepatic function	Retest at a later time (days to weeks) if liver dysfunction is still suspected and patient is anicteric	Consistent with decreased hepatic function
	12–25 µmol/L	Retest at a later time (days to weeks) if liver dysfunction is still suspected and patient is anicteric ¹	Retest at a later time (days to weeks) if liver dysfunction is still suspected and patient is anicteric	Consistent with decreased hepatic function
	>25 µmol/L	Decreased hepatic function possible Retest at a later time if liver dysfunction is still suspected and patient is anicteric ¹	Retest at a later time (days to weeks) if liver dysfunction is still suspected and patient is anicteric	Consistent with decreased hepatic function
		Inadequate fast possible		

¹Preprandial values will occasionally exceed postprandial values. This is usually attributed to spontaneous gall bladder contraction or preprandial specimen obtained prior to a sufficient fast. When preprandial values exceed the postprandial values, always check that pre- and postprandial samples are properly identified.

IDEXX Customer Support
US/Canada/Latin America: 1-800-248-2483
Europe: 0800 1234 3399
New Zealand: 0800-102-084
Australia: 1300 44 33 99
Asia: 0800-291-018

If possible, please have your kit lot number available.

Symbol Descriptions

	Use by date		Authorized Representative in the European Community
	Batch Code (Lot)		Consult instructions for use
	Catalog Number		In vitro diagnostic
	Temperature limitation		
	Manufacturer		

¹SNAP, SNAPshot Dx and VetTest are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2015 IDEXX Laboratories, Inc. All rights reserved • 06-11249-04

IDEXX Laboratories, Inc. One IDEXX Drive, Westbrook, Maine 04092 USA
ECIREP IDEXX Europe B.V. P.O. Box 1334, NL-2130 EK Hoofddorp

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^[*] © 2015 IDEXX Laboratories, Inc. Tutti i diritti sono riservati.